



AlloMap[®] Variability: Long-Term Surveillance for Heart Transplant Patients

 **AlloMap[®]**
Variability Report

Assessing Risk of Future Clinical Events in Heart Transplant Recipients

Need for Personalized Surveillance

Despite progress in the care and management of transplant patients, there remain challenges in optimizing long-term outcomes¹.

Service Offering to Assess Long-Term Risks

CareDx now provides AlloMap Score Variability (AMV) results for patients receiving AlloMap. AMV is a service offering from CareDx that aids in estimating the risk of one or more of the following future clinical events in heart transplant recipients:

- Graft dysfunction²
- Death from any cause
- Re-transplantation

Calculation of AMV

AMV is calculated using the standard deviation of four consecutive, AlloMap Scores. Qualifying AMV scores meet the following criteria:

- ≥ 315 days post transplant
- Time interval between the 1st and 4th score is ≥ 85 and ≤ 780 days.

Clinical Development of AMV

AMV was clinically developed over the course of more than five years and has been assessed in more than 600 patients:

- **Observation** — AlloMap score variability among clinically stable patients without rejection was lower than score variability in patients who experienced rejection; ISHLT³.
- **Hypothesis** — retrospective analysis of randomized trial cohort demonstrated that AMV was predictive of future clinical events; IMAGE⁴.
- **Performance** — evaluation of AMV in an independent patient population demonstrated that variability scores may be predictive of long-term risk; CARGO II⁵.

¹Starling et al. Multicenter Analysis of Immune Biomarkers and Heart Transplant Outcomes: Results of the Clinical Trials in Organ Transplantation-05 Study. *Am J Transplant* 2015 Aug 10 [Epub Ahead of print].

²Graft dysfunction is defined as a clinically significant hemodynamic compromise (left ventricular ejection fraction $\leq 40\%$, use of hemodynamic support such as inotropic medications and/or mechanical assist device, hospitalization for treatment of graft failure or dysfunction).

³Deng, et al. Low Variability of Intraindividual Longitudinal Leukocyte Gene Expression Profiling Cardiac Allograft Rejection Scores. *Transplantation*. 2010; 90(4): 459–461.

⁴Deng, et al. Utility of Gene Expression Profiling Score Variability to Predict Clinical Events in Heart Transplant Recipients. *Transplantation*. 2014; 97(6): 708–714.

⁵Crespo-Leiro, et al. Performance of gene-expression profiling test score variability to predict future clinical events in heart transplant recipients. *BMC Cardiovasc Disord* (2015).



AlloMap Variability Clinical Performance Characteristics⁵

	Variability Score	NPV (95% CI)	% Scores Below	PPV (95% CI)
Lower risk scores	0.1	—	0.0	17.0 (—)
	0.2	100.0 (—)	1.1	17.2 (11.6–23.8)
	0.3	100.0 (—)	3.3	17.8 (12.3–24.7)
	0.4	100.0 (—)	7.7	19.0 (13.1–26.9)
	0.5	100.0 (—)	16.5	22.0 (15.2–30.8)
	0.6	97.0 (91.4–100.0)	24.2	23.3 (15.7–33.0)
	0.7	92.8 (85.4–97.8)	33.0	23.2 (15.4–33.7)
	0.8	91.4 (84.2–96.5)	42.9	25.2 (16.3–37.1)
	0.9	90.3 (83.9–95.2)	56.0	29.7 (18.5–46.0)
	1.0	87.4 (80.1–92.9)	62.6	26.0 (14.8–41.8)
	1.1	87.0 (80.4–92.4)	68.1	27.8 (15.3–46.5)
	1.2	87.0 (80.3–92.2)	71.4	29.9 (16.3–51.0)
Higher risk scores	1.3	86.0 (79.3–91.3)	76.9	29.4 (14.6–53.4)
	1.4	85.3 (78.7–90.6)	83.5	31.9 (13.5–63.2)
	1.5	84.9 (78.5–90.1)	87.9	35.4 (13.5–75.8)
	1.6	83.9 (77.1–89.4)	92.3	29.4 (5.9–100.0)
	1.7	83.7 (77.1–89.1)	94.5	31.9 (0.0–100.0)
	1.8	83.2 (76.6–88.6)	97.8	23.8 (—)
	1.9	83.4 (77.0–88.7)	98.9	100.0 (—)
	2.0	83.4 (77.0–88.7)	98.9	100.0 (—)
	2.1	83.0 (76.6–88.4)	100.0	— (—)


Use of AMV may enable clinicians to personalize patient surveillance and care > 315 days post-transplant. Current performance characteristics associated with AMV suggest that patients with scores⁵:

- ≤ 0.6 may be at a lower risk of future events
- ≥ 1.6 may be at a higher risk of future events

AMV Patient Report Features

- A** Time Post Transplant
- B** Most Recent AlloMap Score
- C** AlloMap NPV and PPV*
- D** Most Recent AMV Score
- E** AMV NPV and PPV
- F** Criteria for Next AMV Score

*NPV and PPV: negative & positive predictive value



To receive AMV contact your CareDx representative or send an email to CareDx Customer Service (customercare@caredx.com).

Upon opting in to receive AMV, a Transplant Center Summary Report (TCSR) will be sent. The TCSR will include the four most recent AlloMap scores and current AMV score (if available) for all active patients in your practice for the current year and previous three calendar years.

You will also receive individual AMV score reports for each patient who receives a traditional AlloMap. Please note that AMV reports will arrive separate from traditional AlloMap Score reports.



CareDx, Inc.
3260 Bayshore Blvd.
Brisbane, CA 94005

1.888.255.6627
www.caredx.com
www.allomap.com
email: inquiries@caredx.com